

# Endonuclease IV, *E. coli*

Cat. No. E70100

Connect with Epicentre on our blog ([epicentral.blogspot.com](http://epicentral.blogspot.com)),  
Facebook ([facebook.com/EpicentreBio](https://facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

## 1. Introduction

Endonuclease IV from *Escherichia coli* is a 32-kD metalloprotein that aids in the repair of damaged DNA. The enzyme functions both as an apurinic/apyrimidinic nuclease<sup>1</sup> and as a 3'-terminal diesterase.<sup>1-4</sup> The latter activity is important in the repair of DNA strand breaks generated by oxidation (e.g., H<sub>2</sub>O<sub>2</sub>) and ionic radiation.<sup>2,3</sup> In such events, the strand breaks terminate with either a 3' phosphate or a deoxyribose fragment, preventing repair by DNA polymerase I or DNA ligase. Endonuclease IV removes the blocking groups, leaving a free 3'-hydroxyl terminus. Although a metalloenzyme, Endonuclease IV is active in the presence of EDTA provided a suitable substrate is present.<sup>4</sup> In addition, the enzyme does not have detectable associated exonuclease or DNA N-glycosylase activities.<sup>1</sup>

Endonuclease IV from *E. coli* is available in a 100-Unit size at a concentration of 2 U/ $\mu$ l.

## 2. Product Specifications

**Storage Temperature:** Store only at -20°C in a freezer without a defrost cycle.

**Storage Buffer:** Endonuclease IV supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 1.0 mM dithiothreitol, 100 mM NaCl, and 0.1% Triton® X-100.

**Unit Definition:** One unit of Endonuclease IV converts 1  $\mu$ g of a partially depurinated, covalently closed supercoiled plasmid DNA to relaxed closed circular form (Form II) in 30 minutes at 37°C.

**Note:** *Partially depurinated plasmid is generated by incubating supercoiled plasmid DNA in 100 mM NaCl and 2.5 mM potassium acetate (pH 4.8) at 70°C for 20 minutes.*

**Quality Control:** The enzyme is function-tested in a reaction containing 50 mM Tris-acetate (pH 7.5), 50 mM KCl, 1 mM EDTA, 1  $\mu$ g partially depurinated pUC19 DNA, and varying amounts of enzyme for 30 minutes at 37°C.

**Contaminating Activity Assays:** Endonuclease IV is free of detectable exo- and endonuclease and RNase activities except for the inherent endonucleolytic properties of the enzyme.

## 3. Related Products

The following products are also available:

- T4 Endonuclease V
- HK™-UNG Thermolabile Uracil N-Glycosylase

## 4. References

1. Ljungquist, S. (1977) *J. Biol. Chem.* **252**, 2808.
2. Demple, B. *et al.*, (1986) *Proc. Natl. Acad. Sci. USA* **83**, 7731.
3. Levin, J.D. *et al.*, (1988) *J. Biol. Chem.* **263**, 8066.
4. Levin, J.D. *et al.*, (1991) *J. Biol. Chem.* **266**, 22893.

*HK is a trademark of Epicentre, Madison, Wisconsin.*

*Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.*

Visit our technical blog: [epicentral.blogspot.com](http://epicentral.blogspot.com)