

EZ-Tn5™ <DHFR-1>Tnp Transposome™ Kit

Cat. No. TSM99D1

Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook (facebook.com/EpicentreBio), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

The EZ-Tn5™ <DHFR-1>Tnp Transposome™ is the stable complex formed between the EZ-Tn5 Transposase enzyme and the EZ-Tn5 <DHFR-1> Transposon. The EZ-Tn5 <DHFR-1> Transposon contains the dihydrofolate reductase (DHFR) gene (trimethoprim resistance) that is functional in *E. coli*, flanked by hyperactive 19-basepair Mosaic End (ME) EZ-Tn5 Transposase recognition sequences. The EZ-Tn5 Transposome can be electroporated into living cells where the EZ-Tn5 Transposase is activated by Mg²⁺ in the host's cellular environment resulting in random insertion of the EZ-Tn5 Transposon into the host genomic DNA.¹⁻⁴

Unlabeled forward and reverse transposon-specific primers are supplied in the kit. The primers can be used for bidirectional DNA sequencing or mapping of transposon insertion sites in target genomic DNAs, including direct sequencing from genomic DNA without cloning.³

2. Applications

- Create gene “knockouts” by insertional mutagenesis, *in vivo*.[†]
- Sequence bacterial genomic DNA directly from transposon insertions, without cloning.
- Insert priming sites for mapping or PCR into chromosomal DNA, *in vivo*.
- Insert a trimethoprim-resistance marker into bacteria

3. Product Specifications

Storage: Store only at –20°C in a freezer without a defrost cycle.

Storage Buffer: The EZ-Tn5 <DHFR-1>Tnp Transposome is supplied in a 50% glycerol solution containing 27.5 mM Tris-HCl (pH 7.5), 50 mM NaCl, 0.3 mM EDTA, 0.05% Triton® X-100, and 0.5 mM dithiothreitol. The DHFR-1 FP-1 Forward and DHFR-1 RP-1 Reverse Primers are supplied in 10 mM Tris-HCl, (pH 7.5), 1 mM EDTA.

Size: Reagents included in the kit are sufficient for 10 *in vivo* transposon insertion reactions.

Quality Control: EZ-Tn5 <DHFR-1>Tnp Transposome activity is assayed by electroporation into a *recA*[–] *E. coli* host strain having a transformation efficiency of >10⁹ cfu/μg DNA. Assays must yield >10⁶ trimethoprim^R colonies/μg or >10⁴ trimethoprim^R colonies/μl of transposome respectively. Primers are function-tested in a DNA cycle sequencing reaction using the SequiTherm EXCEL™ II DNA Sequencing Kit and in a PCR reaction using a plasmid containing an EZ-Tn5 <DHFR-1> Transposon as template.

Contaminating Activity Assays: All components of the EZ-Tn5 <DHFR-1>Tnp Transposome Kit are free of detectable DNase and RNase activities as judged by agarose gel electrophoresis following over-digestion assays, with the exception of the inherent endonucleolytic function of the EZ-Tn5 Transposase.

4. Kit Contents

Desc.	Concentration	Quantity
EZ-Tn5™ <DHFR-1>Tnp Transposome™	@ 15 ng/μl	10 μl
DHFR-1 FP-1 Forward Primer	@ 50 μM	20 μl
DHFR-1 RP-1 Reverse Primer	@ 50 μM	20 μl
Sterile Water		1 ml

5. Related Products

The following products are also available:

- TransforMax™ EC100™ Electrocompetent *E. coli*
- MasterPure™ DNA Purification Kits
- EZ-Tn5™ Tnp Transposome™ Kits
- EZ-Tn5™ Transposase
- EZ-Tn5™ Transposons
- EZ-Tn5™ Insertion Kits

6. Protocol

Electroporation of Host Cells with EZ-Tn5 <DHFR-1>Tnp Transposome and Selection of Transposition Clones:

Electroporate electrocompetent cells using 1 ml of the EZ-Tn5 <DHFR-1>Tnp Transposome. The electrocompetent cells should have a transformation efficiency of $>10^7$ cfu/μg of DNA, but use cells of the highest transformation efficiency possible to maximize the number of transposon insertion clones. Perform electroporation according to the equipment manufacturer's recommendations.


Immediately recover the electroporated cells after electroporation. Even slight delays in initiating the cell recovery process will result in a reduced number of transposition clones. For *E. coli*, add SOC medium to the electroporation cuvette to 1 ml final volume **immediately** after electroporation. Pipette the medium/cells gently to mix. Transfer to a tube and incubate on a 37°C shaker for 30-60 minutes to facilitate cell outgrowth.

Important Note: *EZ-Tn5 <DHFR-1> transposition clones are selected on trimethoprim plates. Trimethoprim is an antimetabolite, not an antibiotic, and the sensitivity of dihydrofolate reductase to this compound varies widely among strains of bacteria, including E. coli. We have found E. coli strain DH10B (Life Technologies, Inc.) to be sensitive to 10 μg/ml of trimethoprim. We recommend testing the sensitivity of your strain by plating on media containing a range (10-200 μg/ml) of trimethoprim. Choose a concentration of that results in no visible growth or a light haze of small, pin point colonies. Another strain must be chosen if growth is confluent.*

If working with *E. coli*, dilute aliquots of the recovered cells (e.g., 1:10 and 1:100). Plate 100 μl of undiluted cells and each cell dilution separately on plates containing trimethoprim (a stock of 10 mg/ml can be made in dimethylformamide). Store the unused portion of the electroporated cells at +4°C for up to 2 days in the event that additional plates need to be prepared. The number of trimethoprim^R colonies/μl of EZ-Tn5 <DHFR-1>Tnp Transposome will be dependent on the transformation efficiency

of the cells used and the level of expression of the DHFR gene in that species. Epicentre's TransforMax EC100 Electrocompetent *E. coli* (available separately) have a transformation efficiency of $>1 \times 10^9$ cfu/ μ g and are ideal for this application.

EZ-Tn5™ <DHFR-1> Transposon
(887 bp.)



Note: Not all restriction enzymes that cut only once are indicated above.

See the following pages for further information.

Primers are not drawn to scale.

DHFR-1 FP-1 Forward Primer	5' GGCGGAAACATTGGATGCGG 3'
DHFR-1 RP-1 Reverse Primer	5' GACACTCTGTTATTACAAATCG 3'
ME = Mosaic End	5' AGATGTGTATAAGAGACAG 3'

Figure 1. EZ-Tn5 <DHFR-1> Transposon.

Restriction Enzymes that cut the EZ-Tn5 <DHFR-1> Transposon one to three times:

Enzyme	Sites	Location	Enzyme	Sites	Location
Acc65 I	1	46	Hinc II	3	189, 679, 839
Acc I	1	838	Hind III	1	855
Aci I	2	732, 750	Hinf I	3	352, 697, 835
Afl III	1	394	HinP I	1	734
AlwN I	1	297	Hpa I	2	189, 679
Apo I	1	34	Hpa II	2	51, 571
Ase I	1	132	Hpy188 I	2	196, 414
BamH I	1	825	HpyCH4 III	1	141
Ban I	1	46	HpyCH4 IV	2	396, 425
Ban II	1	44	HpyCH4 V	3	66, 845, 853
Bfa I	1	832	Kpn I	1	50
BfuA I	2	790, 850	Mae II	2	396, 425
Bpu10 I	1	795	Mbo II	3	212, 249, 790
BsaH I	1	769	Mly I	1	844
BsaJ I	1	271	Msp I	2	52, 572
BseY I	1	601	MspA1 I	2	152, 752
BsiE I	1	389	Mwo I	2	721, 742
BsiHKA I	1	44	Nci I	1	51
Bsl I	1	799	Nco I	1	271
BsmA I	2	8, 875	Nde I	1	774
Bsp1286 I	1	44	Nla IV	2	48, 827
BspD I	2	29, 527	Nsp I	3	66, 172, 853
BspM I	2	790, 850	Pac I	1	136
Bsr I	1	616	Ple I	1	843
BssK I	3	49, 736, 757	PspG I	2	736, 757
BstDS I	1	275	Pst I	1	847
BstF5 I	1	823	Rsa I	1	48
BstN I	2	738, 759	Sac I	1	44
BstU I	1	734	Sal I	1	837
BstY I	1	825	Sau96 I	1	259
Btg I	1	271	Sbf I	1	847
Cla I	2	29, 527	ScrF I	3	51, 738, 759
Dde I	2	195, 795	SfaN I	3	155, 441, 808
Dra I	2	100, 303	Sfc I	1	843
Dsa I	1	271	Sph I	3	66, 172, 853
EcoR I	1	34	Sty I	1	271
Fau I	2	145, 814	Tfi I	2	352, 697
Fnu4H I	2	689, 732	Tse I	1	688
Hae I	1	615	Tsp45 I	2	155, 709
Hae III	2	260, 615	Tsp4C I	1	141
Hha I	1	736	Xba I	1	831

Restriction Enzymes that cut the EZ-Tn5 <DHFR-1> Transposon four or more times:

Alu I	CviJ I	Mae III	Mse I	Taq I
Alw I	Dpn I	Mbo I	Nla III	Tsp509 I
Cac8 I	Hph I	Mnl I	Sau3A I	

Restriction Enzymes that do not cut the EZ-Tn5 <DHFR-1> Transposon:

Aat II	Bme1580 I	BstZ17 I	NgoM IV	Sfi I
Acl I	BmgB I	Bsu36 I	Nhe I	Sfo I
Afe I	Bmr I	Bts I	Not I	SgrA I
Afl II	Bsa I	Dra III	Nru I	Sim I
Age I	BsaA I	Drd I	Nsi I	Sma I
Ahd I	BsaB I	Eae I	PaeR7 I	Sml I
Ale I	BsaW I	Eag I	Pci I	SnaB I
Apa I	BsiW I	Ear I	PfiF I	Spe I
ApaB I	Bsm I	Eco47 III	PfiM I	Srf I
ApaL I	BsmB I	EcoN I	Pme I	Sse8647 I
Asc I	BspE I	EcoO109 I	Pml I	Ssp I
AsiS I	BspH I	EcoR V	PpuM I	Stu I
Ava I	BspLU11 I	Fse I	PshA I	Swa I
Ava II	BsrB I	Fsp I	Psi I	Tat I
Avr II	BsrD I	Gdi II	PspOM I	Tli I
Bbs I	BsrF I	Hae II	Pvu I	TspR I
BbvC I	BsrG I	Hpy99 I	Pvu II	Tth111 I
BciV I	BssH II	Mfe I	Rsr II	Xcm I
Bcl I	BssS I	Mlu I	Sac II	Xho I
BfrB I	BstAP I	Msc I	SanD I	Xma I
Bgl I	BstB I	Msl I	Sap I	Xmn I
Bgl II	BstE II	Nae I	Sca I	
Blp I	BstX I	Nar I	SexA I	

EZ-Tn5™ <DHFR-1> Transposon 887 bp.

1	CTGTCTCTTA	TACACATCTC	AACCATCATC	GATGAATTCG	AGCTCGGTAC
51	CCGGATAGAC	GGCATGCACG	ATTTGTAATA	ACAGAGTGTC	TTGTATTTTT
101	AAAGAAAGTC	TATTTAATAC	AAGTGATTAT	ATTAATTAAC	GGTAAGCATC
151	AGCGGGTGAC	AAAACGAGCA	TGCTTACTAA	TAAAATGTTA	ACCTCTGAGG
201	AAGAATTGTG	AAACTATCAC	TAATGGTAGC	TATATCGAAG	AATGGAGTTA
251	TCGGGAATGG	CCCTGATATT	CCATGGAGTG	CCAAAGGTGA	ACAGCTCCTG
301	TTTAAAGCTA	TTACCTATAA	CCAATGGCTG	TTGGTTGGAC	GCAAGACTTT
351	TGAATCAATG	GGAGCATTAC	CCAACCGAAA	GTATGCGGTC	GTAACACGTT
401	CAAGTTTTAC	ATCTGACAAT	GAGGACGTAT	TGATCTTTCC	ATCAATTAAA
451	GATGCTTTAA	CCAACCTAAA	GAAAATAACG	GATCATGTCA	TTGTTCAGG
501	TGGTGGGGAG	ATATACAAAA	GCCTGATCGA	TCAAGTAGAT	ACACTACATA
551	TATCTACAAT	AGACATCGAG	CCGGAAGGTG	ATGTTTACTT	TCCTGAAATC
601	CCCAGCAATT	TTAGGCCAGT	TTTTACCCAA	GACTTCGCCT	CTAACATAAA
651	TTATAGTTAC	CAAATCTGGC	AAAAGGGTTA	ACAAGTGGCA	GCAACGGATT
701	CGCAAACCTG	TCACGCCTTT	TGTGCCAAAA	GCCGCGCCAG	GTTTGCGATC
751	CGCTGTGCCA	GGCGTTAGGC	GTCATATGAA	GATTTCCGGT	ATCCCTGAGC
801	AGGTGGCGGA	AACATTGGAT	GCGGGGATCC	TCTAGAGTCG	ACCTGCAGGC
851	ATGCAAGCTT	CAGGGTTGAG	ATGTGTATAA	GAGACAG	

The transposon sequence can be downloaded at www.epicentre.com/sequences

Primer Information

DHFR-1 FP-1 Forward Primer

5' - GGCGGAAACATTGGATGCGG - 3'

Length: 20 nucleotides

G+C content: 12

Molecular Weight: 6,228 daltons

Temperatures of Dissociation & Melting:

T_d : 75°C (nearest neighbor method)

T_m : 72°C (% G+C method)

T_m : 64°C ([2 (A+T) + 4 (G+C)] method)

T_m : 64°C ((81.5 + 16.6 (log [Na⁺])) + ([41 (#G+C) - 500] / length) method)
where [Na⁺] = 0.1 M

DHFR-1 RP-1 Reverse Primer

5' - GACTACTGTATTACAAATCG - 3'

Length: 22 nucleotides

G+C content: 8

Molecular Weight: 6,672 daltons

Temperatures of Dissociation & Melting:

T_d : 56°C (nearest neighbor method)

T_m : 66°C (% G+C method)

T_m : 60°C ([2 (A+T) + 4 (G+C)] method)

T_m : 57°C ((81.5 + 16.6 (log [Na⁺])) + ([41 (#G+C) - 500] / length) method)
where [Na⁺] = 0.1 M

7. References

1. Hoffman, L.M. *et al.*, (1999) *Current Genetics* **35**, 304.
2. Reznikoff, W. and Goryshin, I.Y. (1999) *Epicentre Forum* **6** (2), 5.
3. Hoffman, L.M. and Jendrisak, J. (1999) *Epicentre Forum* **6** (3), 1.
4. Goryshin, I.Y. *et al.*, (2000) *Nat. Biotechnol.* **18**, 97.

EC100, EZ-Tn5, HyperMu, MasterPure, SequiTherm EXCEL, Transposome, and TransforMax are trademarks of Epicentre, Madison, Wisconsin. Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

¹Use of Transposome™ Complexes for in vivo Transposon Insertion, including, but not limited to EZ-Tn5™ and HyperMu™ Transposome™ complexes, is covered by U.S. Patent Nos. 6,159,736 and 6,294,385; European Patent No. 1115856 and other patents issued or pending, exclusively licensed to Epicentre. These products are accompanied by a limited non-exclusive license for nonprofit organizations to use the purchased products solely for life science research. Specifically, researchers at nonprofit organizations may use these products in their research and they may transfer derivatives of the products to colleagues at other nonprofit organizations provided that such colleagues agree in writing to be bound by the terms and conditions of this label license. Researchers may not transfer these products or its derivatives to researchers at other organizations that are not nonprofit organizations without the express written consent of Epicentre and without those entities having an express license from Epicentre for the use of the products. Other than as specifically provided here, an authorized transferee of these products shall have no right to transfer these products or its derivatives to any other person or entity. A nonprofit organization performing research using this product for a for-profit organization shall also be considered for-profit. For-profit entities purchasing these products shall have three (3) months to evaluate the products in their research. Thereafter, such for-profit entities shall either contact Epicentre and take a license for their continued use of the products or destroy the products and all derivatives thereof. Please contact Epicentre with respect to licenses for commercial use, including manufacturing, therapeutic or diagnostic use.

Visit our technical blog: epicentral.blogspot.com