

EZ-Tn5™ <TET-1> Insertion Kit

Cat. No. EZI921T

Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook (facebook.com/EpicentreBio), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

Transposons are mobile DNA sequences found in the genomes of prokaryotes and eukaryotes. Transposon tagging has long been recognized as a powerful research tool for randomly distributing primer binding sites, creating gene “knockouts”, and introducing a physical tag or a genetic tag into large target DNAs. One frequently used transposition system is the Tn5 system isolated from gram-negative bacteria.

Though a naturally occurring transposition system, the Tn5 system can be readily adapted for routine use in research laboratories for the following reasons:

- 1) Tn5 transposase is a small, single subunit enzyme that has been cloned and purified to high specific activity.
- 2) Tn5 transposase carries out transposition without the need for host cell factors.
- 3) Tn5 transposon insertions into target DNA are highly random.
- 4) Tn5 transposition proceeds by a simple “cut and paste” process. Although the chemistry is unique, the result is similar to using a restriction endonuclease, with random sequence specificity, accompanied by a DNA ligase activity.
- 5) Tn5 transposase will transpose **any** DNA sequence contained between its short 19 basepair Mosaic End (ME) Tn5 transposase recognition sequences.

In 1998 Goryshin and Reznikoff¹ demonstrated that a fully functional Tn5 transposition system could be reconstituted *in vitro*. Additionally, the transposition efficiency of this system has been increased more than 1,000-fold compared to wild-type Tn5 by introducing mutations in the transposase gene and in the 19-bp Tn5 ME transposase recognition sequence.

Epicentre’s EZ-Tn5 Transposon Tools (kits and reagents) are based on the hyperactive Tn5 transposition system developed by Goryshin and Reznikoff.

2. Kit Contents

Desc.	Concentration	Quantity
Reagents included in the kit are sufficient for 10 <i>in vitro</i> transposon insertion reactions.		
EZ-Tn5™ Transposase in 50% glycerol containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, 0.1 mM EDTA, 1 mM dithiothreitol, and 0.1% Triton® X-100.	10 U @ 1 U/μl	10 μl
EZ-Tn5™ <TET-1> Transposon: 1 pmol in TE Buffer (10 mM Tris-HCl [pH 7.5] and 1 mM EDTA).	@ 0.1 pmol/μl (0.11 μg/μl)	10 μl
EZ-Tn5™ 10X Reaction Buffer: 0.50 M Tris-acetate (pH 7.5), 1.5 M potassium acetate, 100 mM magnesium acetate, and 40 mM spermidine.		100 μl
EZ-Tn5™ 10X Stop Solution: 1% SDS solution.		100 μl
TET-1 FP-1 Forward Primer: 1 nmol in TE Buffer (10 mM Tris-HCl [pH 7.5] and 1 mM EDTA).	@ 50 μM	20 μl
TET-1 RP-1 Reverse Primer: 1 nmol in TE Buffer (10 mM Tris-HCl [pH 7.5] and 1 mM EDTA).	@ 50 μM	20 μl
pUC19/3.4 Control Target DNA: 1 μg A 3.4-kb <i>Hpa</i> II fragment of bacteriophage DNA cloned into the <i>Acc</i> I site of pUC19, in TE Buffer (10 mM Tris-HCl [pH 7.5] and 1 mM EDTA).	@ 0.1 μg/μl	10 μl
Sterile Water:		1 μl

3. Related Products

The following products are also available:

- TransforMax™ EC100™ Electrocompetent *E. coli*
 - Transformation efficiency >10⁹ cfu/μg DNA.
- Colony Fast-Screen™ Kit
- MasterPure™ DNA Purification Kits
- EZ-Tn5™ Tnp Transposome™ Kits
- EZ-Tn5™ Transposase
- EZ-Tn5™ Transposons
- EZ-Tn5™ Insertion Kits

How the EZ-Tn5 <TET-1> Insertion Kit Works.

The EZ-Tn5 <TET-1> Insertion Kit can be used to randomly insert primer binding sites and a Tetracycline resistance selection marker into target DNA *in vitro*. A single 2-hour *in vitro* reaction randomly inserts the <TET-1> Transposon into the target DNA. Use an aliquot of the reaction to transform *E. coli* such as Epicentre's TransforMax™ EC100™ Electrocompetent *E. coli* strain and select on Tetracycline plates. Only those clones harboring DNA containing the <TET-1> Transposon will grow.

The EZ-Tn5 <TET-1> Insertion Kit can be used for:

- Faster sequencing of large DNA molecules, as compared to primer walking, random subcloning, or generating nested deletions with exonuclease III and mung bean nuclease.
- Making insertion mutants or gene “knockouts” *in vitro*.
- Introducing a Tetracycline resistance selection marker into any DNA.

Fig. 1 describes the steps involved when using the EZ-Tn5 <TET-1> Insertion Kit. The process can be summarized as follows:

Preparation

- Prepare 0.2 μg of recombinant DNA for the EZ-Tn5 <TET-1> insertion reaction.

Day 1

- Perform the 2-hour *in vitro* EZ-Tn5 <TET-1> insertion reaction.
- Transform competent *recA*⁻ *E. coli* with 1 μl of the reaction mix.
- Select for Tetracycline-resistant transposon insertion clones on Tetracycline plates overnight.

Day 2

- Prepare DNA from Tetracycline-resistant colonies.
- (Optional) Map the EZ-Tn5 <TET-1> Transposon insertion sites.
- (Optional) DNA sequence chosen clones bidirectionally using the unlabeled forward and reverse transposon-specific primers supplied in the kit.




Figure 1. EZ-Tn5 <TET-1> Transposon Insertion Protocol.

4. Materials

Storage

Store EZ-Tn5 Insertion Kits only at -20°C in a freezer without a defrost cycle. The enzyme solution will not freeze. Other component solutions will freeze. After thawing, be sure to mix the contents of each tube thoroughly before using. After use, return all of the kit components to -20°C for storage.

Performance Specifications and Quality Control

The EZ-Tn5 <TET-1> Insertion Kit is function-tested by performance of the kit's *in vitro* control reaction followed by electroporation into a *recA*⁻ *E. coli* host strain having a transformation efficiency of $>10^9$ cfu/ μg DNA. Transposition frequency, defined as the ratio of the number of Tet^R clones divided by the number of transformants resistant to the antibiotic marker of the target vector, (Tet^R colonies/Amp^R colonies; for the control DNA) must be $>0.5\%$ (commonly at 10%) and transposition efficiency must be $>10^6$ Tet^R colonies/ μg target DNA. Primers are function-tested in a DNA cycle sequencing reaction using the SequiTherm EXCEL™ II DNA Sequencing Kit and in a PCR reaction using a plasmid containing an EZ-Tn5 <TET-1> Transposon as template. All components of the EZ-Tn5 <TET-1> Insertion Kit are free of detectable DNase and RNase activities as judged by agarose gel electrophoresis following over-digestion assays, with the exception of the inherent endonucleolytic function of the EZ-Tn5 Transposase.

5. Transposon Insertion Reaction

Target DNA Preparation

The target DNA **must not** contain a Tetracycline resistance gene. The transposon insertion reaction is not significantly inhibited by high levels of RNA contamination in target DNA preparations. However, if the target DNA is heavily contaminated with chromosomal DNA, which is a direct competitor for target transposition, the number of clones will be greatly reduced.

Plasmid and cosmid clones can be purified by standard minilysate procedures and used as target DNA in the insertion reaction without further clean-up. Low copy-number vectors, for example BAC or fosmid clones, are often contaminated with a higher molar proportion of *E. coli* chromosomal DNA thus reducing the transposon insertion frequency. Therefore, BAC and fosmid DNA should be purified, to remove the chromosomal DNA prior to the insertion reaction.

In Vitro Transposon Insertion Reaction

Reaction conditions have been optimized to maximize the efficiency of the EZ-Tn5 Transposon insertion while minimizing multiple insertion events. **Be sure to calculate the moles of target DNA used in the reaction and add an equimolar amount of the EZ-Tn5 <TET-1> Transposon.** If necessary, dilute the EZ-Tn5 <TET-1> Transposon with TE buffer (10 mM Tris-HCl [pH 7.5], 1 mM EDTA).

1. Prepare the transposon insertion reaction mixture by adding in the following order:

1	µl	EZ-Tn5 10X Reaction Buffer
0.2	µg	target DNA*
x	µl	molar equivalent EZ-Tn5 <TET-1> Transposon
x	µl	sterile water to a reaction volume of 9 µl
1	µl	EZ-Tn5 Transposase
<hr/>		
10	µl	Total reaction volume
2. Incubate the reaction mixture for 2-hours at 37°C.
3. Stop the reaction by adding 1 µl EZ-Tn5 10X Stop Solution. Mix and heat for 10 minutes at 70°C. Proceed to section 4, **Transformation and Recovery**, or store the reaction mixture at -20°C.

* Calculation of µmol target DNA:

$$\mu\text{mol target DNA} = \mu\text{g target DNA} / [(\# \text{ base pairs in target DNA}) \times 660]$$

For example: 0.2 µg of control pUC19/3.4 DNA which is 6,100 bp

$$= 0.2 \mu\text{g} / [6,100 \text{ bp} \times 660] = 0.05 \times 10^{-6} \mu\text{mol} = 0.05 \text{ pmol}$$

6. Selection of Transposon Insertion Clones

Transformation and Recovery

The number of EZ-Tn5 Transposon insertion clones obtained per reaction depends on, among other factors, the transformation efficiency of the competent cells used. The greater the transformation efficiency of the competent cells, the greater the number of insertion clones obtained. We recommend using electrocompetent or chemically competent *recA*⁻ *E. coli* with a transformation efficiency of >10⁸ cfu/µg of DNA. A *recA*⁻ strain of *E. coli* is recommended to eliminate the possibility of generating multimeric forms of the vector. Finally, the host strain used **must not** have a tetracycline resistance marker when used with the EZ-Tn5 <TET-1> Transposon. Epicentre's TransforMax EC100 Electrocompetent *E. coli* (available separately) have a transformation efficiency of >1 x 10⁹ cfu/µg and are ideal for this application.

- 1) Using 1 µl of the insertion reaction mixture, transform *recA*⁻ *E. coli*. If electrocompetent cells are used, perform electroporation according to the equipment manufacturer's recommended conditions. Use of >1 µl of the transposon insertion reaction for transformation may result in arcing. The unused portion of the transposon insertion reaction can be stored at -20°C for future use.
- 2) Recover the electroporated cells by adding medium (**Important: medium must not contain Mg²⁺**) to the electroporation cuvette to 1 ml final volume immediately after electroporation. Pipette the medium/cells gently to mix. Transfer to a tube and incubate on a 37°C shaker for 30-60 minutes to facilitate cell outgrowth.

Plating and Selecting Transformants

EZ-Tn5 <TET-1> Transposon insertion clones are selected on tetracycline-containing plates.

Important: Plate medium must not contain Mg^{2+} .

Epicentre recommends storing tetracycline-containing plates in the dark, and using plates less than one month old.

- 1) If transformation was done using cells with an efficiency of $>5 \times 10^8$ cfu/ μ g DNA, it may be necessary to dilute the cells 1:10 or 1:100 prior to plating. Plate portions (e.g., 100 μ l) of cells onto LB plates containing 10 μ g/ml tetracycline. Store the unused portion of the recovered cells at +4°C for up to 2 days in the event additional plates need to be prepared.
- 2) (Optional) To determine the transposon insertion efficiency, plate identical dilutions and dilution aliquots of the transformation reaction on a second plate containing an antibiotic specific for selecting target DNA (e.g., 100 μ g/ml ampicillin for the control DNA). The transposition frequency is given by the ratio of Tet^R/Amp^R clones for the control DNA.
- 3) Grow plates overnight at 37°C. Assuming a transposon insertion efficiency of 1% and use of high purity target DNA (i.e., little or no chromosomal DNA contamination), one should see 100-500 Tet^R clones per plate. If too few (or too many) colonies appear, replat the transformed cells at a lower (or higher) dilution.

Typical results obtained with the EZ-Tn5 <TET-1> Insertion Kit are:

- Transposon insertion frequency = 0.5-20%
- Transposition clones per μ g target DNA = 1×10^5 - 1×10^8
- Transposition clones per 10 μ l transposon insertion reaction = 1×10^4 - 1×10^7

The actual number of EZ-Tn5 <TET-1> insertion clones obtained will vary depending on factors such as target DNA size and the transformation efficiency of the competent cells used to recover the transposon insertion clones.

7. DNA Sequencing of Transposon Insertion Clones

Transposon Insertion Mapping (optional)

EZ-Tn5 Transposon insertion clones can be sequenced bidirectionally using the unlabeled forward and reverse transposon-specific primers provided in the kit. The insertion site of each clone can also be mapped prior to sequencing, depending on the cost-effectiveness and experience of the user.

EZ-Tn5 <TET-1> Transposon insertion sites can be mapped by size analysis of PCR products using colony minilysate DNA as a template. To map the insertion sites, use the TET-1 FP-1 or TET-1 RP-1 primers provided with the kit and a vector-specific flanking primers (not provided).

Alternatively, insertion sites can be mapped by restriction endonuclease digest(s). Use the nucleotide sequence and restriction information of the EZ-Tn5 <TET-1> Transposon provided in the Appendix for reference.

Primer Considerations

The TET-1 FP-1 Forward and TET-1 RP-1 Reverse Primers supplied with the kit have been constructed to minimize homology to commonly used cloning vectors. **However, the sequence of each primer should be compared to that of the user's specific cloning vector to ensure minimal sequence homology to the vector.** The sequence and theoretical melting temperatures for each primer are presented in the Appendix.

Note: The sequence of the TET-1 FP-1 Forward and TET-1 RP-1 Reverse Primers do not function well as IRD800-labeled sequencing primers. We recommend using the TET-1 FP-3 Forward and TET-1 RP-3 Reverse Primer sequences for this purpose. See the Appendix.

Note: Occasionally a clone will yield the sequence of the cloning vector. This occurs when the EZ-Tn5 Transposon randomly inserts into a non-essential region of the vector rather than into the DNA insert. The frequency of this occurrence is dependent on the size of the DNA insert relative to the size of non-essential regions of vector. The larger the DNA insert, the less frequently an insertion will occur into the vector.

Target Site Duplication

EZ-Tn5 Transposase-catalyzed transposon insertion results in the generation of a 9-bp target site sequence duplication where one copy immediately flanks each side of the inserted transposon. This is important to consider when assembling the nucleotide sequence of a recombinant clone insert. The process of transposon insertion site duplication is depicted in Fig. 2.

Distinguishing Transposon Sequence from Insert Sequence

Since the primers provided in the EZ-Tn5 <TET-1> Insertion Kit anneal to a region near the ends of the transposon, the first sequence data obtained from each sequencing reaction is that of EZ-Tn5 <TET-1> Transposon DNA. The sequence of the 19-bp EZ-Tn5 Transposase recognition sequence (ME) found at the junction of the transposon and the target clone insert DNA (present in all insertion clones), is a useful landmark to distinguish vector sequence from target clone insert sequence (see also Fig. 3).




Figure 2. EZ-Tn5 Transposase Insertion Site Duplication Process.

- 2) **The sequencing primer used has significant homology to the cloning vector or to the DNA being sequenced.** Check homology of primer against the vector. Alter primer annealing conditions or synthesize a new primer with less homology.
- 3) **Components of the DNA sequencing kit and/or of the electrophoresis step are compromised.** Reverify the integrity of the components of the kit and/or electrophoresis step with appropriate controls. Use a new kit and/or new reagents.

Confluent plates following transformation

- 1) **Target DNA or host cells have the same selective marker (antibiotic resistance) as the trans-poson used.** Use a different host for transformation and retransform with a portion of your remaining reaction.

Some host cells carry Tet^R transposons. Confirm that the genotype of the host strain chosen for the transformation is not Tet^R.

9. Appendix

Primer Data

TET-1 FP-1 Forward Primer

5' - GGGTGCGCATGATCCTCTAGAGT - 3'

Length: 23 nucleotides

G+C content: 13

Molecular Weight: 7,074 daltons

Temperatures of Dissociation & Melting:

T_d: 71°C (nearest neighbor method)

T_m: 75°C (% G+C method)

T_m: 72°C ([2 (A+T) + 4 (G+C)] method)

T_m: 66°C ((81.5 + 16.6 (log [Na⁺])) + ([41 (#G+C) - 500] / length) method)
where [Na⁺] = 0.1 M

TET-1 RP-1 Reverse Primer

5' - TAAATTGCACTGAAATCTAGAAATA - 3'

Length: 25 nucleotides

G+C content: 6

Molecular Weight: 7,657 daltons

Temperatures of Dissociation & Melting:

T_d: 61°C (nearest neighbor method)

T_m: 64°C (% G+C method)

T_m: 62°C ([2 (A+T) + 4 (G+C)] method)

T_m: 55°C ((81.5 + 16.6 (log [Na⁺])) + ([41 (#G+C) - 500] / length) method)
where [Na⁺] = 0.1 M

Note: Use the following alternative primer sequences for IRD800-labeled sequencing primers.

TET-1 FP-3 Forward Primer

5' - GCATCTCGGGCACGTTGGGTCCT - 3'

Length: 23 nucleotides

G+C content: 15

Molecular Weight: 7,905 daltons

Temperatures of Dissociation & Melting:

T_d : 80°C (nearest neighbor method)

T_m : 79°C (% G+C method)

T_m : 76°C ([2 (A+T) + 4 (G+C)] method)

T_m : 70°C $((81.5 + 16.6 (\log [\text{Na+}])) + ([41 (\#G+C) - 500] / \text{length}) \text{ method})$
where $[\text{Na+}] = 0.1 \text{ M}$

TET-1 RP-3 Reverse Primer

5' - GGGGCTGACTTCAGGTGCTACATT - 3'

Length: 25 nucleotides

G+C content: 13

Molecular Weight: 8,561 daltons

Temperatures of Dissociation & Melting:

T_d : 73°C (nearest neighbor method)

T_m : 76°C (% G+C method)

T_m : 76°C ([2 (A+T) + 4 (G+C)] method)

T_m : 66°C $((81.5 + 16.6 (\log [\text{Na+}])) + ([41 (\#G+C) - 500] / \text{length}) \text{ method})$
where $[\text{Na+}] = 0.1 \text{ M}$


EZ-Tn5 <TET-1> Transposon Sequence

EZ-Tn5™ <TET-1> Transposon 1,674 bp.

1	CTGTCTCTTA	TACACATCTC	AACCATCATC	GATGAATTCG	AGCTCGGTAC	CCGGGGATCA
61	TCTTATTAAT	CAGATAAAAT	ATTTCTAGAT	TTCAGTGCAA	TTTATCTCTT	CAAATGTAGC
121	ACCTGAAGTC	AGCCCCATAC	GATATAAGTT	GTAATTCTCA	TGTTTGACAG	CTTATCATCG
181	ATAAGCTTTA	ATGCGGTAGT	TTATCACAGT	TAAATTGCTA	ACGCAGTCAG	GCACCGTGTA
241	TGAAATCTAA	CAATGCGCTC	ATCGTCATCC	TCGGCACCGT	CACCCTGGAT	GCTGTAGGCA
301	TAGGCTTGGT	TATGCCGGTA	CTGCCGGGCC	TCTTGCGGGA	TATCGTCCAT	TCCGACAGCA
361	TCGCCAGTCA	CTATGGCGTG	CTGCTAGCGC	TATATGCGTT	GATGCAATTT	CTATGCGCAC
421	CCGTTCTCGG	AGCACTGTCC	GACCCTTTG	GCCGCCGCC	AGTCTGTCTC	GCTTCGCTAC
481	TTGAGGCCAC	TATCGACTAC	GCGATCATGG	CGACCACACC	CGTCTGTGG	ATCCTCTACG
541	CCGACGCAT	CGTGGCCGGC	ATCACCCGGC	CCACAGGTGC	GGTTGTGTCG	GCCTATATCG
601	CCGACATCAC	CGATGGGGAA	GATCGGGCTC	GCCACTTCGG	GCTCATGAGC	GCTTGTTCG
661	GCGTGGGTAT	GGTGGCAGGC	CCCGTGGCCG	GGGACTGTT	GGCGCCATC	TCCCTGCATG
721	CACCATTCCCT	TGCGGCGGCG	GTGTCACAACG	GCCTCAACCT	ACTACTGGGC	TGCTTCCTAA
781	TGCAGGAGTC	GCATAAGGGA	GAGCGTCGAC	CGATGCCCTT	GAGAGCCTTC	AACCCAGTCA
841	GCTCCTCCG	GTGGGCGCGG	GGCATGACTA	TCGTCGCCGC	ACTTATGACT	GTCCTCTTTA
901	TCATGCAACT	CGTAGGACAG	GTGCCGGCAG	CGCTCTGGGT	CATTTTCGGC	GAGGACCGCT
961	TTCGCTGGAG	CGCGACGATG	ATCGGCCTGT	CGCTTGCGGT	ATTCCGGAATC	TTGCACGCC
1021	TCGCTCAAGC	CTTCGTCACT	GGTCCCGCCA	CCAAACGTTT	CGCGAGAAG	CAGGCCATTA
1081	TCGCCGCAT	GGCGGCCGAC	GCGCTGGGCT	ACGTCTTGCT	GGCGTTCGCG	ACGCGAGGCT
1141	GGATGGCCTT	CCCCATTATG	ATTCTTCTCG	CTTCCGGCGG	CATCGGGATG	CCCCTGTGTC
1201	AGGCCATGCT	GTCCAGGCAG	GTAGATGACG	ACCATCAGGG	ACAGCTTCAA	GGATCGCTCG
1261	CGGCTCTTAC	CAGCCTAACT	TCGATCATTG	GACCGCTGAT	CGTCACGGCG	ATTTATGCCG
1321	CCTCGGCAG	CACATGGAAC	GGGTTGGCAT	GGATTGTAGG	CGCCGCCCTA	TACCTTGCT
1381	GCCTCCCCGC	GTTGCGTTCGC	GGTGATGGA	GCCGGGCCAC	CTCGACCTGA	ATGGAAGCCG
1441	GCGCACCTC	GCTAACGGAT	TCACCCTCC	AAGAATTGGA	GCCAATCAAT	TCTTGCGGAG
1501	AACTGTGAAT	GCGCAAACCA	ACCCTTGGA	GAACATATCC	ATCGCGTCCG	CCATCTCCAG
1561	CAGCGCACGC	GGCGCATCTC	GGGCACGTTG	GGTCTTGCC	ACGGGTGCGC	ATGATCTCT
1621	AGAGTCGACC	TGCAGGCATG	CAAGCTTCAG	GGTTGAGATG	TGTATAAGAG	ACAG

The transposon sequence can be downloaded at the URL: <http://www.epicentre.com/sequences>

EZ-Tn5 <TET-1> Transposon Schematic Map

EZ-Tn5™ <TET-1> Transposon
(1,674 bp)

Note: Not all restriction enzymes that cut only once are indicated above.

See the following pages for further information.

Hind III cuts twice in the map above.

Primers are not drawn to scale.

TET-1 FP-1 Forward Primer

5' GGGTGC GCATGATCCTCTAGAGT 3'

TET-1 RP-1 Reverse Primer

5' TAAATTGCACTGAAATCTAGAAATA 3'

ME = Mosaic End

5' AGATGTGTATAAGAGACAG 3'

Figure 4. EZ-Tn5 <TET-1> Transposon

EZ-Tn5 <TET-1> Transposon Restriction Data**Restriction Enzymes that cut the EZ-Tn5 <TET-1> Transposon one to three times:**

Enzyme	Sites	Location	Enzyme	Sites	Location
Acc65 I	1	46	EcoR V	1	341
Acc I	2	806, 1625	Fsp I	3	416, 1512, 1608
Acl I	1	1055	Hinc II	2	807, 1626
Afe I	3	388, 650, 931	Hind III	2	183, 1642
ApaB I	1	1206	HpyCH4 IV	3	1055, 1111, 1585
Apo I	1	34	Kpn I	1	50
Ase I	1	66	Mae II	3	1055, 1111, 1585
Ava I	2	50, 1578	Mly I	2	795, 1631
BamH I	1	529	Msc I	1	1598
Ban II	3	44, 629, 643	Mse I	3	66, 188, 210
Bbs I	1	884	Msl I	1	1185
Bfa I	3	85, 384, 1619	MspA1 I	1	1295
BfuA I	2	1208, 1637	Nhe I	1	383
Bgl I	2	1089, 1323	Nru I	1	1128
Bme1580 I	1	1585	Nsp I	2	720, 1640
Bmr I	3	453, 774, 828	PflM I	2	1475, 1524
BsaW I	1	847	Ple I	2	794, 1630
BseY I	1	1103	PpuM I	1	1591
BsiE I	3	443, 810, 1096	PshA I	1	870
Bsm I	1	1513	PspG I	3	283, 1212, 1593
BsmA I	2	8, 1662	Pst I	1	1634
BspD I	2	29, 178	Rsa I	2	48, 319
BspH I	1	643	Sac I	1	44
BspM I	2	1208, 1637	Sal I	2	805, 1624
BstAP I	1	1205	Sbf I	1	1634
BstDS I	2	686, 1603	Sfc I	2	292, 1630
BstN I	3	285, 1214, 1595	SgrA I	1	564
BstY I	1	529	Sim I	2	938, 1591
Btg I	2	682, 1599	Sma I	1	52
Cla I	2	29, 178	Sml I	2	818, 1024
Dsa I	2	682, 1599	Sph I	2	720, 1640
Eag I	1	1093	Ssp I	1	80
Ear I	1	112	Sty I	1	1523
Eco47 III	3	388, 650, 931	Tfi I	3	1006, 1160, 1458
EcoN I	1	780	TspR I	3	99, 444, 1048
EcoO109 I	2	678, 1591	Xba I	2	84, 1618
EcoR I	1	34	Xma I	1	50

Restriction Enzymes that cut the EZ-Tn5 <TET-1> Transposon four or more times:

Aci I	BssK I	Hae III	Mbo I	Sau3A I
Alu I	BstF5 I	Hha I	Mbo II	Sau96 I
Alw I	BstU I	Hinf I	Mnl I	ScrF I
Ava II	Cac8 I	HinP I	Msp I	SfaN I
Ban I	CviJ I	Hpa II	Mwo I	Sfo I
BsaH I	Dpn I	Hph I	Nae I	Taq I
BsaJ I	Eae I	Hpy188 I	Nar I	Tse I
BsiHKA I	Fau I	Hpy99 I	Nci I	Tsp45 I
Bsl I	Fnu4H I	HpyCH4 III	NgoM IV	Tsp4C I
Bsp1286 I	Gdi II	HpyCH4 V	Nla III	Tsp509 I
Bsr I	Hae I	Mae III	Nla IV	
BsrF I	Hae II			

Restriction Enzymes that do not cut the EZ-Tn5 <TET-1> Transposon:

Aat II	Bgl II	BstB I	Not I	Sca I
Afl II	Blp I	BstE II	Nsi I	SexA I
Afl III	BmgB I	BstX I	Pac I	Sfi I
Age I	Bpu10 I	BstZ17 I	Paer7 I	SnaB I
Ahd I	Bsa I	Bsu36 I	Pci I	Spe I
Ale I	BsaA I	Bts I	PfIF I	Srf I
AlwN I	BsaB I	Dde I	Pme I	Sse8647 I
Apa I	BsiW I	Dra I	Pml I	Stu I
ApaL I	BsmB I	Dra III	Psi I	Swa I
Asc I	BspE I	Drd I	PspOM I	Tat I
AsiS I	BspLU11 I	Fse I	Pvu I	Tli I
Avr II	BsrB I	Hpa I	Pvu II	Tth111 I
BbvC I	BsrD I	Mfe I	Rsr II	Xcm I
BciV I	BsrG I	Mlu I	Sac II	Xho I
Bcl I	BssH II	Nco I	SanD I	Xmn I
BfrB I	BssS I	Nde I	Sap I	

10. References

Cited:

1. Goryshin, I.Y. and Reznikoff, W.S. (1998) *J. Biol. Chem.* **273**, 7367.

Related:

1. Goryshin, I.Y. and Reznikoff, W.S. (1998) *J. Biol. Chem.* **273**, 7367.
2. Zhou, M. *et al.*, (1998) *J. Mol. Biol.* **276**, 913.
3. Zhou, M. and Reznikoff, W.S. (1997) *J. Mol. Biol.* **271**, 362.
4. Mahnke-Braam, I.A. and Reznikoff, W.S. (1998) *J. Biol. Chem.* **273**, 10908.

These products are intended for research use only.

EZ-Tn5™ Transposon Tools for in vitro transposon insertion are covered by U.S. Patent Nos. 5,925,545; 5,948,622; 5,965,443, and 6,437,109; European Patent No. 0927258, and other patents issued or pending, exclusively licensed or assigned to Epicentre. These products are accompanied by a limited non-exclusive license for the purchaser to use the purchased product(s) solely for in vitro transposon insertion for life science research. Purchase of these products does not grant rights to: (1) offer products, components of products, or any derivatives thereof for resale; or (2) to distribute or transfer the products, components of products, or any derivatives thereof to third parties. Contact Epicentre for information on licenses for uses other than life science research

EZ-Tn5, TransforMax, EC100, SequiTherm EXCEL, Fast-Screen, MasterPure, and Transposome are trademarks of Epicentre, Madison, Wisconsin.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

Visit our technical blog: epicentral.blogspot.com