

# MaxPlax™ Lambda Packaging Extracts

Cat. Nos. MP5105, MP5110, and MP5120

Connect with Epicentre on our blog ([epicentral.blogspot.com](http://epicentral.blogspot.com)),  
Facebook ([facebook.com/EpicentreBio](https://facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

## 1. Introduction

MaxPlax™ Lambda Packaging Extracts are a convenient, high-efficiency system designed for *in vitro* lambda packaging reactions. MaxPlax Lambda Packaging Extracts are supplied as predispensed single-tube reactions that have been optimized for packaging of methylated and unmethylated DNA. The packaging extracts routinely yield packaging efficiencies of  $> 1 \times 10^9$  pfu/ $\mu$ g of Control  $\lambda$  DNA. The extracts can be used in the construction of representative cDNA libraries and genomic cloning of highly modified (methylated) DNA into  $\lambda$ -phage or cosmid vectors.

Traditional packaging extracts are derived from two complementary lysogenic *E. coli* strains, BHB2690 and BHB2688, as described by Hohn (1979).<sup>1</sup> The MaxPlax extracts utilize a new packaging strain, NM759\*, reported by Gunther, Murray and Glazer (1993).<sup>2</sup> This strain, which replaces strain BHB2690 in the preparation of the sonication extract, is a restriction-free K12-derived strain deficient in the production of  $\lambda$ -phage capsid protein D. When combined with the complementary freeze-thaw extract from strain BHB2688\*\*,<sup>1</sup> deficient in the production of  $\lambda$ -phage capsid protein E, an extremely high-efficiency of packaging for  $\lambda$  DNA is obtained. Moreover, the ability to package  $\lambda$  DNA bearing the mammalian methylation pattern is greatly enhanced, as evidenced by the high efficiency of  $\lambda$ -vector rescue from transgenic mouse DNA.<sup>2</sup> The lack of restriction activity has been shown to be crucial for the high efficiency rescue of lambda shuttle vectors from transgenic mouse DNA.<sup>2,3</sup>

\*NM759: [W3110 *recA56*,  $\Delta$ (*mcrA*) *e14*,  $\Delta$ (*mrr-hsd-mcr*), (*aimm434*, *clts*, *b2*, *red3*, *Dam15*, *Sam7*)/ $\lambda$ ]

\*\*BHB2688: [N205 *recA*\*, (*aimm434* *clts*, *b2*, *red3*, *Eam4*, *Sam7*)/ $\lambda$ ]

## 2. Kit Contents

Desc.	Concentration	Quantity
MaxPlax Lambda Packaging Extracts are available in 5-10-, and 20-extract sizes.		
The 20-extract size kit contains the following reagents:		
MaxPlax Lambda Packaging Extracts	50 $\mu$ l/tube	20 tubes
Ligated Lambda Control DNA ( $\lambda$ c1857 <i>Sam7</i> )	(1 $\mu$ g @ 0.02 $\mu$ g/ $\mu$ l)	50 $\mu$ l
Control plating <i>E. coli</i> strain LE392MP, glycerol stock		250 $\mu$ l
[F <sup>-</sup> <i>e14</i> (McrA <sup>-</sup> ) $\Delta$ ( <i>mcrC-mrr</i> ) (Tet <sup>R</sup> ) <i>hsdR514 supE44 supF58 lacY1</i> or $\Delta$ ( <i>lacIZY</i> )6 <i>galK2 galT22 metB1 trpR55</i> $\lambda$ -]		

**Note:** MaxPlax Lambda Packaging Extracts are supplied as freeze-thaw/sonicate extracts in unlabeled single tubes. The extracts, Ligated Lambda Control DNA, and Control Plating Strain LE392MP are packaged together in a CO<sub>2</sub>-impermeable foil pouch.

**Store the MaxPlax Lambda Packaging Extracts at -70°C or below. Exposure to higher temperature will decrease packaging efficiencies.**

### 3. Product Specifications

**Storage:** Store the control host bacteria and the MaxPlax Lambda Packaging Extracts at  $-70^{\circ}\text{C}$ . Exposure to higher temperatures will greatly compromise packaging extract efficiency. Avoid long term storage of product in the presence of dry ice. Once removed from the foil package, avoid any exposure to dry ice. Store the remainder of the kit components at  $-20^{\circ}\text{C}$ . After thawing, store the Control DNA at  $4^{\circ}\text{C}$ .

**Storage Buffers:** MaxPlax Lambda Packaging Extracts are supplied as unlabeled single tubes of freeze-thaw/sonicate extracts. Control plating bacteria are supplied as a glycerol stock. Control ligated  $\lambda$  DNA is supplied in 1X Ligation Buffer.

**Guaranteed Stability:** MaxPlax Lambda Packaging Extracts are guaranteed to maintain a packaging efficiency of  $>1.0 \times 10^9$  pfu/ $\mu\text{g}$  of control  $\lambda$  DNA, when stored as directed for one year from the date of purchase.

### 4. Protocol for Packaging Lambda DNA

This protocol can be used for the positive control reaction as well as for experimental reactions. The positive control reactions must be plated on the control host bacterial strain (LE392MP) included with the MaxPlax Extracts. The proper bacterial plating strain for the experimental reactions will vary depending on the cloning vector used. See the vector manufacturer's recommendations for the proper strain and plating media requirements. Ligation reactions may be added directly to the packaging extracts. When doing so, it is important to: a) add a volume of 10  $\mu\text{l}$  or less to the packaging reaction, and b) heat inactivate the ligase (i.e., treatment at  $65^{\circ}\text{C}$  for 15 minutes) as active DNA ligase will decrease packaging efficiencies.

#### Solutions:

Phage Dilution Buffer	LB Broth (1 Liter)	LB Plates
10 mM Tris-HCl (pH 8.3)	10 g Bacto-tryptone	LB Broth with 1.5% (w/v)
100 mM NaCl	5 g Bacto-yeast extract	Bacto-agar
10 mM $\text{MgCl}_2$	10 g NaCl	<b>LB Top Agar</b>
	Adjust pH to 7.0 with NaOH	LB Broth with 0.7% (w/v)
		Bacto-agar

#### Plating Bacteria Preparation:

1. The day before performing the packaging reactions, inoculate 50 ml of supplemented (10 mM  $\text{MgSO}_4$ ) LB broth with a single colony of the plating bacterial strain and shake overnight at  $37^{\circ}\text{C}$ .
2. The day of the packaging reactions, inoculate 50 ml of supplemented (10 mM  $\text{MgSO}_4$  + 0.2% maltose) LB broth with 5 ml of the overnight culture and shake at  $37^{\circ}\text{C}$  to an  $\text{OD}_{600} = 0.8-1.0$ . Store the cells at  $4^{\circ}\text{C}$  until needed; cells may be stored for up to 72 hours.

### Packaging Reactions:

1. Thaw the appropriate number of packaging extracts at room temperature. For every two packaging reactions, thaw one extract then place on ice.
2. When thawed, immediately transfer half (25  $\mu$ l) of each packaging extract to a second 1.5-ml tube and place on ice.
3. Add the substrate DNA (10  $\mu$ l [0.2  $\mu$ g] of the control DNA) to a tube containing 25  $\mu$ l of extract. If performing an odd number of packaging reactions, the remaining 25  $\mu$ l of extract can be refrozen at  $-70^{\circ}\text{C}$ .
4. Mix by pipetting several times; avoid the introduction of air bubbles. Return all of the contents to the bottom of the tube by brief centrifugation if necessary.
5. Incubate the reaction(s) at  $30^{\circ}\text{C}$  for 90 minutes.
6. At the end of this incubation, add the additional 25  $\mu$ l of thawed extract to each reaction tube at  $30^{\circ}\text{C}$  (If performing two packaging reactions, thaw another tube of extract and add 25  $\mu$ l to each tube.) and incubate the reaction(s) for an additional 90 minutes at  $30^{\circ}\text{C}$ .
7. Add 500  $\mu$ l of phage dilution buffer and mix by gentle vortexing.
8. Add 25  $\mu$ l of chloroform and mix by gentle vortexing (store at  $4^{\circ}\text{C}$ ).
9. Assay the packaged phage by titering on the appropriate bacterial strain (LE392MP for the control).

### Titering Phage Extracts:

1. Make serial dilutions of the packaged phage in phage dilution buffer. Use  $10^{-5}$  and  $10^{-6}$  dilutions for the control reactions.  
10<sup>-2</sup> dilution is 10  $\mu$ l of packaged phage particles into 990  $\mu$ l of phage dilution buffer; vortex mix.  
10<sup>-4</sup> dilution is 10  $\mu$ l of 10<sup>-2</sup> dilution into 990  $\mu$ l phage dilution buffer; vortex mix.  
10<sup>-5</sup> dilution is 100  $\mu$ l of 10<sup>-4</sup> dilution into 900  $\mu$ l phage dilution buffer; vortex mix.  
10<sup>-6</sup> dilution is 10  $\mu$ l of 10<sup>-4</sup> dilution into 990  $\mu$ l phage dilution buffer; vortex mix.
2. Add 100  $\mu$ l of the appropriate serial dilutions to 100  $\mu$ l of prepared plating bacteria (use LE392MP for the control reactions) and incubate for 15 minutes at  $37^{\circ}\text{C}$ .
3. Add 3.0 ml of melted supplemented (10 mM  $\text{MgSO}_4$ ) LB top agar (cooled to  $\sim 48^{\circ}\text{C}$ ). Vortex gently and pour onto pre-warmed ( $37^{\circ}\text{C}$ ) LB plates. Allow the top agar to solidify and then incubate overnight at  $37^{\circ}\text{C}$ .
4. Count the plaques and determine the titer (pfu/ml) and packaging efficiency (See sample calculations).

**Sample Calculations:**

If there were 110 plaques on a  $10^{-6}$  dilution plate, then the titer, pfu/ml, (where pfu represents plaque forming units) of this reaction would be:

$$\frac{(\# \text{ of plaques}) (\text{dilution factor}) (1000 \mu\text{l/ml})}{(\text{volume of phage plated } [\mu\text{l}])} \text{ OR } \frac{(110 \text{ pfu}) (10^6) (1000 \mu\text{l/ml})}{(100 \mu\text{l})} = 1.1 \times 10^9 \text{ pfu/ml}$$

The packaging efficiency (pfu/ $\mu\text{g}$  DNA) of this reaction would be:

$$\frac{(\# \text{ of plaques}) (\text{dilution factor}) (\text{total reaction vol.})}{(\text{vol. of dilution plated}) (\text{amount of DNA packaged})} \text{ OR } \frac{(110 \text{ pfu}) (10^6) (550 \mu\text{l})}{(100 \mu\text{l}) (0.2 \mu\text{g})} = 3 \times 10^9 \text{ pfu}/\mu\text{g}$$

**5. Related Products**

The following products are also available:

Desc.	Concentration	Quantity
<b>pWEB-TNC™ Cosmid Cloning Kit</b> WEBC931		1 Kit
<b>pWEB™ Cosmid Cloning Kit</b> PC8805		1 Kit
<b>EpiFOS™ Fosmid Library Production Kit</b> FOS0901		1 Kit
<b>Lambda Terminase</b> LT4450	2 U/ $\mu\text{l}$	50 Units
LT44200	2 U/ $\mu\text{l}$	200 Units

**6. References**

1. Hohn, E.G. (1979) *Methods Enzymol.* **68**, 299.
2. Gunther, E.G. *et al.*, (1993) *Nucl. Acids Res.* **21**, 3903.
3. Kohler, S.W. *et al.*, (1990) *Nucl. Acids Res.* **18**, 3007.

*EpiFOS*, *MaxPlax*, *pWEB*, and *pWEB-TNC* are trademarks of Epicentre, Madison, Wisconsin.

Visit our technical blog: [epicentral.blogspot.com](http://epicentral.blogspot.com)