

Rec J Exonuclease, *E. coli*

Cat. Nos. RJ411050 and RJ411250

Important! *The storage temperature for Rec J Exonuclease has changed.
The new storage temperature is –65°C to –85°C.*

Epicentre (an Illumina company) warrants that its products will retain full activity for 1 year from the date of receipt by the user if used and stored properly. For full warranty details, see Terms and Conditions available on the Epicentre website at www.epibio.com.

Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook (facebook.com/EpicentreBio), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

Rec J Exonuclease, derived from *E. coli*, catalyzes removal of deoxyribonucleotide monophosphates from single stranded DNA (ssDNA) in a 5'→3' direction. The enzyme requires magnesium (Mg²⁺) and can be heat inactivated by incubation at 65°C for 20 minutes.

The enzyme is supplied with a 10X Reaction Buffer.

Applications

- Removal of primers from completed PCR reactions.
- Degradation of linear ssDNA in double-stranded DNA (dsDNA) and plasmid preps.

2. Product Specifications

Storage: Store only at –65°C to –85°C in a freezer without a defrost cycle.

Storage Buffer: Rec J Exonuclease is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1.0 mM dithiothreitol (DTT), 0.1 mM EDTA and 0.1% Triton® X-100.

Unit Definition: One unit of Rec J Exonuclease results in the acid-solubilization of 1 nmol of nucleotides from activated single-stranded calf thymus DNA in 30 minutes at 37°C.

10X Rec J Exonuclease Reaction Buffer: 330 mM Tris-acetate (pH 7.5), 660 mM potassium acetate, 100 mM magnesium acetate and 5 mM DTT.

Quality Control: Rec J Exonuclease is function-tested in a reaction containing 33 mM Tris-acetate (pH 7.5), 66 mM potassium acetate, 10 mM magnesium acetate, 0.5 mM DTT, 10 µg of heat-denatured activated calf thymus DNA and varying amounts of Rec J Exonuclease.

Contaminating Activity Assays: Rec J Exonuclease is free of detectable RNase, DNA endonuclease, and dsDNA exonuclease activities.

3. Related Products

The following products are also available:

- Exonuclease I, *E. coli*
- Exonuclease III, *E. coli*
- Exonuclease VII, *E. coli*
- Lambda Exonuclease
- Rec BCD Nuclease
- T5 Exonuclease
- Terminator™ 5'-Phosphate-Dependent Exonuclease
- Mung Bean Nuclease
- OmniCleave™ Endonuclease

OmniCleave and Terminator are trademarks of Epicentre, Madison, Wisconsin.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

Visit our technical blog: epicentral.blogspot.com