

## RNase-Free DNase I

Cat. Nos. D9902K, D9905K, and D9910K

Connect with Epicentre on our blog ([epicentral.blogspot.com](http://epicentral.blogspot.com)),  
Facebook ([facebook.com/EpicentreBio](https://facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

## 1. Introduction

RNase-Free DNase I is an endonuclease that efficiently hydrolyzes double- (ds) or single-stranded DNA to a mixture of short oligo- and mononucleotides. In the presence of  $Mg^{2+}$ , cleavage of each strand of a dsDNA substrate proceeds independently.<sup>1</sup> In contrast, in the presence of  $Mn^{2+}$ , the enzyme cleaves both strands of DNA at approximately the same site to generate molecules with blunt ends or 1- or 2-base overhangs<sup>1</sup> that can be blunted with T4 DNA Polymerase.

Provided at a concentration of 1 MBU/ $\mu$ l, RNase-Free DNase I is available in 2,500-, 5,000-, and 10,000-MBU sizes and includes 10X Reaction Buffer.

## 2. Applications

- Elimination of the DNA template following *in vitro* RNA synthesis with T7 Phage RNA Polymerase.
- Characterization of DNA:protein interactions by of "DNase I footprinting".<sup>1,2</sup>
- Treatment of RNA prior to RT-PCR.<sup>3</sup>
- Radiolabeling of DNA by nick translation.<sup>1,4</sup>

## 3. Product Specifications

**Storage:** Store only at  $-20^{\circ}C$  in a freezer without a defrost cycle.

**Storage Buffer:** RNase-Free DNase I is supplied in a 50% glycerol solution containing 10 mM Tris-HCl (pH 7.5), 50 mM  $CaCl_2$ , and 10 mM  $MgCl_2$ .

**DNase I 10X Reaction Buffer:** 100 mM Tris HCl (pH 7.5), 25 mM  $MgCl_2$ , and 5 mM  $CaCl_2$ .

**Unit Definition:** One Molecular Biology Unit (MBU) of RNase-Free DNase I produces an increase in the  $A_{260}$  of a solution of dsDNA, of 0.001 per minute at  $25^{\circ}C$ . Functionally, 1 MBU completely digests 1  $\mu$ g of pUC19 DNA to oligonucleotides in 10 minutes at  $37^{\circ}C$ .

**Quality Control:** RNase-Free DNase I is function-tested in two assay systems. A hyperchromicity assay is performed in a reaction containing 50  $\mu$ g/ml native calf thymus DNA, 0.1 M sodium acetate (pH 5.0), 5 mM  $MgCl_2$ , 2 mM  $CaCl_2$ , and varying amounts of enzyme. A digestion assay is performed in a reaction containing 33 mM Tris-acetate (pH 7.5), 66 mM potassium acetate, 10 mM magnesium acetate, 0.5 mM dithiothreitol, 2 mM  $CaCl_2$ , 1.0  $\mu$ g of pUC19 DNA, and varying amounts of enzyme.

**Contaminating Activity Assays:** RNase-Free DNase I is free of detectable RNase activities as assayed by PAGE analysis of 1  $\mu$ g of a synthetic RNA transcript following incubation with enough RNase-Free DNase I to completely digest 1 mg of DNA.

## 4. Example Protocol

1. Dilute DNase I 10X Reaction Buffer to 1X using RNase-Free water.
2. Prepare 50 µl of a working DNase I Solution for each sample to be treated by adding 5 µl of RNase-Free DNase I to 45 µl of 1X Reaction Buffer (from Step 1).
3. Completely resuspend 5 µg of a nucleic acid pellet in 50 µl of working DNase I solution.
4. Incubate at 37°C for 10 minutes.

## 5. References

1. Sambrook, J. *et al.*, (1989) in: *Molecular Cloning: A Laboratory Manual (2nd ed.)*, Cold Spring Harbor Laboratory Press, New York.
2. Galas, D.J. and Schmitz, A. (1978) *Nucleic Acids Res.* **5**, 3157.
3. Kienzle, N. *et al.*, (1996) *BioTechniques* **20**, 612.
4. Rigby, P.W.J. *et al.*, (1977) *J. Mol. Biol.* **113**, 237.

## 6. Related Products

Cat. #	Concentration	Quantity
<b>Plasmid-Safe™ ATP-Dependent DNase</b>		
E3101K	10 U/µl	1,000 Units
E3105K	10 U/µl	5,000 Units
E3110K	10 U/µl	10,000 Units
<b>Exonuclease I, <i>E. coli</i></b>		
X40501K	20 U/µl	1,000 Units
X40505K	20 U/µl	5,000 Units
X40520K	20 U/µl	20,000 Units
<b>Exonuclease III, <i>E. coli</i></b>		
EX4405K	200 U/µl	5,000 Units
EX4425K	200 U/µl	25,000 Units
<b>Mung Bean Nuclease</b>		
M8202K	50 U/µl	2,000 Units
M8205K	50 U/µl	5,000 Units
<b>OmniCleave™ Endonuclease</b>		
OC7810K	200 U/µl	10,000 Units
OC7850K	200 U/µl	50,000 Units

*OmniCleave* and *Plasmid-Safe* are trademarks of Epicentre, Madison, Wisconsin.

Visit our technical blog: [epicentral.blogspot.com](http://epicentral.blogspot.com)