

T7 RNA Polymerase

For all sizes of T7 RNA Polymerase
(Cat. Nos. T7905K, T7950K, TM910K, TH950K, TU950K)

Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook ([facebook.com/EpicentreBio](https://www.facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

In vitro transcription with T7 RNA Polymerase enables the synthesis of specific RNA from DNA sequences cloned downstream of the appropriate promoter. RNA can be radiolabeled for use as a probe in various blotting and *in situ* hybridization techniques. Nonradioactive RNA synthesized with T7 RNA Polymerase can be used in studies of translation, RNA processing, and antisense RNA.

2. Product Specifications

Storage: Store only at -20°C in a freezer without a defrost cycle.

Storage Buffer: T7 RNA Polymerase is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, 0.1 mM EDTA, 1 mM DTT, and 0.1% Triton[®] X-100.

Unit Definition: One unit converts 1 nmol of ribonucleoside triphosphates (NTPs) into acid-insoluble material in 60 minutes at 37°C .

5X Transcription Buffer: 200 mM Tris-HCl (pH 7.5), 30 mM MgCl_2 , 50 mM NaCl, and 10 mM spermidine. DTT and NTPs must also be added to the final reaction.

Activity Assay: The unit assay is performed in a reaction containing 40 mM Tris-HCl (pH 7.5), 6 mM MgCl_2 , 10 mM NaCl, 2 mM spermidine, 10 mM DTT, 5 μg DNA template, 0.5 mM of each NTP, and varying amounts of enzyme.

Contaminating Activity Assays: T7 RNA Polymerase is free of detectable DNA exonuclease and endonuclease, RNase, and *E. coli* RNA polymerase activities.

3. Example Protocol

3.A. Standard Transcription Reaction: Synthesis of Nonradioactive RNA

- Combine the following reaction components at room temperature in the order given:

	<u>Final Concentration</u>
x μl RNase-Free water	---
0.5-1 μg linearized template DNA with appropriate promoter	25-50 ng/ μl
4 μl 5X Transcription Buffer	1X
1 μl 10 mM ATP	0.5 mM
1 μl 10 mM CTP	0.5 mM
1 μl 10 mM GTP	0.5 mM
1 μl 10 mM UTP	0.5 mM
2 μl 100 mM DTT	10 mM
10 U T7 RNA Polymerase	0.5 U/ μl
<hr/>	
20 μl Total reaction volume	

- Incubate at 37°C for 2 hours.
- (Optional) Add 1 μl (1 MBU) of RNase-Free DNase I and incubate for 15 minutes at 37°C .

The volume of this reaction may be increased or decreased proportionately.

3.B. Synthesis of Radioactive RNA

Radioactive RNA may be synthesized by adding a radioactive nucleoside triphosphate and decreasing the final concentration of the corresponding nonradioactive NTP to 15 μM . As the total concentration of the labeled NTP will be limiting, virtually all of the radioactivity will be incorporated into RNA. For example, to label with α - ^{32}P -CTP, first prepare a stock of 100 μM CTP by diluting the 10 mM CTP stock 1 to 100 with RNase-Free H_2O . Then add 3 μl of the 100 μM CTP and the desired amount of α - ^{32}P -CTP in place of the 10 mM CTP and some of the water in the above protocol.

4. Related Products

The following products are also available:

Cat. #	Concentration	Quantity
AmpliScribe™ T7-Flash™ Transcription Kit		
ASF3257		25 Reactions
ASF3507		50 Reactions
DuraScribe® T7 Transcription Kit		
DS010910		10 Reactions
DS010925		25 Reactions
RNase-Free DNase I		
D9902K	1 U/ μl	2,500 MBU
D9905K	1 U/ μl	5,000 MBU
D9910K	1 U/ μl	10,000 MBU (2 x 5,000 MBU)

DuraScribe is a registered trademark of, and AmpliScribe, and T7-Flash are trademarks of Epicentre, Madison, Wisconsin.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

Visit our technical blog: epicentral.blogspot.com