

# Tobacco Acid Pyrophosphatase

Cat. Nos. T19050, T19100, T19250, T19500, and T81050

Discontinued

## 1. Introduction

Tobacco Acid Pyrophosphatase (TAP) hydrolyzes various pyrophosphate bonds, including those in adenosine triphosphate (ATP), cyclic nucleotides, and dinucleotides, but not those in RNA or DNA.<sup>1</sup> In addition, the enzyme cleaves the pyrophosphate bond of the 5'-terminal methylated guanine nucleotide "cap" of eukaryotic messenger RNAs.<sup>2</sup> The cap structure is also found on many small nuclear RNAs (snRNAs), heterogeneous nuclear RNAs, and many viral RNAs. Complete hydrolysis results in removal of the  $\beta$  and  $\gamma$  phosphates, leaving only the  $\alpha$  phosphate attached. The resulting 5'-monophosphorylated terminus may be ligated to a 3'-hydroxylated terminus using T4 RNA Ligase or dephosphorylated with APex™ Heat-Labile Alkaline Phosphatase for end labeling.

TAP is available in 50-, 100-, 250-, and 500-Unit sizes at a concentration of 10 U/ $\mu$ l. The 50-Unit size is also available at 5 U/ $\mu$ l. A 10X TAP Reaction Buffer is provided with the enzyme.

## 2. Applications

**Cloning of Full-Length cDNAs.** Rapid amplification of cDNA ends (RACE) is a method whereby unknown sequences 5' and 3' of a known sequence are obtained through reverse transcription, followed by amplification and cloning.<sup>3</sup> RNA is decapped and an oligonucleotide anchor sequence is ligated to the 5' end using T4 RNA Ligase.

**Radiolabeling RNA.**<sup>4</sup> Following removal of the cap structure with TAP and dephosphorylation with APex Phosphatase, RNA can be end-labeled with T4 Polynucleotide Kinase and  $\gamma$ -[<sup>32</sup>P]-ATP and used for sequencing or as a hybridization probe.

**RNA 5'-End Conversion.** In addition to decapping capped mRNA, TAP can also be used to convert 5'-triphosphate RNA into 5'-monophosphate RNA.

## 3. Product Specifications

**Storage:** Store only at -20°C in a freezer without a defrost cycle.

**Storage Buffer:** TAP is supplied in a 50% glycerol solution containing 10 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1.0 mM dithiothreitol, 0.1 mM EDTA, and 0.01% Triton® X-100.

**Unit Definition:** One unit releases 1 nmol of inorganic phosphate from m<sup>7</sup>GpppG in 30 minutes at 37°C under standard assay conditions.

**Note:** One Epicentre m<sup>7</sup>GpppG Unit is approximately equivalent to 15 Epicentre ATPase Units.

**10X TAP Reaction Buffer:** 500 mM sodium acetate (pH 6.0), 10 mM EDTA, 1%  $\beta$ -mercaptoethanol, and 0.1% Triton X-100.

**Contaminating Activity Assays:** TAP is free of detectable RNase activity as judged by incubation of 1  $\mu$ g of an SP6 transcript with 50 units of TAP for 1 hour at 37°C, followed by agarose gel electrophoresis. The enzyme is also free of detectable exonuclease and endonuclease activities.

## 4. Suggested Protocols

**For Ligation of Oligonucleotides to RNA:** Perform Part A; proceed to the ligation reaction.

**For End Labeling of RNA:** Perform Parts A, B, and C.

### 4.A. Removal of the Cap from RNA

1. Assemble the reaction in a microcentrifuge tube on ice in the order given (see Note).

x	μl	Nuclease-Free Water
25	pmol	Capped RNA
5	μl	10X TAP Reaction Buffer
25	U	TAP
<hr/>		
50	μl	Total reaction volume

2. Incubate at 37°C for 1-2 hours.

**Note:** TAP does not require ATP for activity. *m*<sup>7</sup>GpppG is likewise not required for enzyme activity, but is used as a substrate in the measurement of enzyme activity. One picomole of a 1-kb RNA is approximately equivalent to 330 ng.

### 4.B. Removal of the α Phosphate from RNA

1. Equilibrate the reaction from Part 4.A at 30°C and assemble the following reaction:

50	μl	TAP-treated RNA reaction mixture from Part 4.A
38	μl	Nuclease-Free Water
10	μl	10X APex Alkaline Phosphatase Buffer
2	μl	APex Alkaline Phosphatase
<hr/>		
100	μl	Total reaction volume

2. Incubate at 37°C for 1 hour.
3. Inactivate the APex Phosphatase by extracting with an equal volume of phenol:chloroform (1:1, v/v). Transfer the aqueous phase to a clean tube. Re-extraction of the organic phase with buffer may improve recovery of the RNA. Extract the combined aqueous phases with an equal volume of chloroform. Ethanol-precipitate and resuspend the RNA in sterile Nuclease-Free Water.

### 4.C. End-Labeling Dephosphorylated RNA with T4 Polynucleotide Kinase (T4 PNK)

1. Assemble the reaction in a microcentrifuge tube on ice in the order given.

25	pmol	dephosphorylated RNA
25	pmol	γ-[ <sup>32</sup> P]-ATP (3,000-7,000 Ci/mmol)
x	μl	10X T4 PNK Buffer
5-10	U	T4 PNK

2. Incubate at 37°C for 30 minutes.
3. Inactivate the T4 PNK by heating the reaction mixture at 65°C for 5 minutes. Precipitate the labeled RNA with ethanol.

## 5. References

1. Shinshi, H. *et al.*, (1976) *Biochemistry* **15**, 2185.
2. Shinshi, H. *et al.*, (1976) *FEBS Letters* **65**, 254.
3. Schaefer, B.C. (1995) *Anal. Biochem.* **227**, 255.
4. Efstratiadis, A. *et al.*, (1977) *Nucleic Acids Res.* **4**, 4165.

## 6. Related Products

The following products are also available:

Cat. #	Concentration	Quantity
<b>APex™ Heat-Labile Alkaline Phosphatase</b> AP49100	1 Reaction/ $\mu$ l	100 Reactions
<b>T4 Polynucleotide Kinase</b> P0503K	10 U/ $\mu$ l	3,000 Units
<b>RNA 5' Polyphosphatase</b> RP8092H	20 U/ $\mu$ l	200 Units

*APex is a trademark of Epicentre, Madison, Wisconsin.*

*Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.*

Discontinued