

Uracil-DNA Excision Mix

Cat. No. UEM04100

Connect with Epicentre on our blog (epicentral.blogspot.com),
Facebook (facebook.com/EpicentreBio), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

1. Introduction

Uracil-DNA Excision Mix* contains two enzymes, HK™-UNG (Heat-Killable Uracil N-Glycosylase [UNG]) and Endonuclease IV. HK-UNG cleaves the uracil base from a uracil-deoxynucleotide in any DNA, creating an abasic site at the location of dUTP incorporation. Endonuclease IV then cleaves the phosphodiester bond at this abasic site generating a series of digestion fragments. The size of the fragments defines the presence and location of the uracil residues in the treated DNA.

Uracil-DNA Excision Mix (100 µl) is provided with 100 µl of a 10X Uracil Excision Enzyme Buffer.

2. Product Specifications

Storage: Store only at –20°C in a freezer without a defrost cycle.

Uracil Excision Enzyme 10X Buffer: is 500 mM Tris-HCl (pH 9.0), 200 mM $[\text{NH}_4]_2\text{SO}_4$ and 100 mM EDTA.

Activity Assay: Uracil Excision Enzyme Mix is function-tested in an assay by incubating DNA containing at least one uracil residue, in a standard reaction and analyzing the digestion products by polyacrylamide gel electrophoresis.

Contaminating Activity Assays: The Uracil Excision Enzyme Mix is free of detectable non-specific DNA exonuclease and endonuclease, and RNase activities.

3. Related Products

The following products are also available:

- MasterPure™ Nucleic Acid Purification Kits
- BuccalAmp™ DNA Extraction Kits
- dNTP Solutions
- dUTP Solution

4. General Considerations for Optimization of Excision/Cleavage Reactions

1. **Amount of DNA:** The optimal amount of DNA to be digested in the excision reaction depends upon the specific activity of the labeled primer and the length and its nucleotide composition. Excision/cleavage reactions containing DNA labeled with fluorescent dyes or nonradioactive molecules require more product for adequate detection (2-3 pmol, see below). Likewise, target sequences composed of a large percentage of A+T will require more DNA in the reaction.
2. **Excision/cleavage product concentration:** If necessary, concentrate the excision/cleavage reaction products by ethanol precipitation before loading the samples on a polyacrylamide gel.

* U.S. Patent No. 6,190,865.

Excision/Cleavage of DNA Containing dUTP

DNA containing uracil residues is cleaved using the Uracil-DNA Excision Mix and then subsequently subjected to denaturing gel electrophoresis.

A. Assembly of the Excision/Cleavage Reaction

The volumes listed below are for a single excision reaction containing a total reaction volume of 10 μ l. This volume may be increased or decreased by proportionally adjusting the individual reaction components.

1. Thaw and thoroughly mix all of the reagents before dispensing; keep on ice.

Combine the following reagents, on ice:

x μ l	deionized water
2.5-5 μ l	DNA containing dUTP (0.75-3 pmol)
1 μ l	10X Uracil Excision Enzyme Buffer
1 μ l	Uracil-DNA Excision Mix
<hr/>	
10 μ l	Total reaction volume

2. Incubate for 30 minutes at 37°C.
3. Add 5 μ l of a Stop/Loading Buffer (ie. 95% formamide, 10 mM EDTA, 0.1% xylene cyanol, 0.1% bromophenol blue [final pH 9.5]) to each reaction; adjust the volume proportionally if the total reaction volume is different from 10 μ l.

Note: Xylene cyanol and bromophenol blue may interfere with some fluorescence detectors. When using such detectors, use a stop/loading buffer containing 95% formamide that is compatible with the detection system (alternative dye).

4. Cleavage products may be stored at -20°C for up to two weeks.

B. Typical Gel Electrophoresis Procedure

1. Prepare a standard 6-8% polyacrylamide sequencing gel containing 8 M urea.
2. Spin the tubes briefly in a microcentrifuge.
3. Heat the tubes 5 minutes at 75°C.
4. Load 1-3 μ l (50-200 fmol of Uracil-DNA Excision Mix product) per well.

BuccalAmp, HK, and MasterPure are trademarks of Epicentre, Madison, Wisconsin.

Visit our technical blog: epicentral.blogspot.com