

Mung Bean Nuclease

Cat. Nos. **M8202K** and **M8205K**

Mung Bean Nuclease (EC 3.1.30.x) is a single-strand-specific nuclease purified from sprouts of mung bean *Vigna radiata*. The enzyme degrades single-stranded DNA or RNA to nucleoside-5'-monophosphates, but does not digest double-stranded DNA, double-stranded RNA, or DNA / RNA hybrids.¹⁻³ Mung Bean Nuclease is preferable to S1 nuclease for most applications because it has lower intrinsic activity on duplex DNA.⁴ The extremely low exonuclease activity of Mung Bean Nuclease is demonstrated in high resolution transcript mapping experiments with end-labeled single-stranded DNA probes.^{3,5}

Applications

- High resolution mapping of termini and exon structure of RNA transcripts ("Berk-Sharp" or "S1-mapping" experiments).^{3,5-8}
- Restriction site modification or removal by digestion of protruding single-stranded ends.³
- Unidirectional deletion of DNA sequences (in combination with Exonuclease III) to facilitate subcloning in kilobase-sequencing strategies (i.e., an alternative to "Bal 31" deletions).⁹
- Removal of hairpin structures during cDNA synthesis.¹⁰

Product Specifications

Storage Temperature: Store only at -20°C in a freezer without a defrost cycle.

Storage Buffer: Mung Bean Nuclease is supplied in a 50% glycerol solution containing 10 mM Tris-HCl (pH 7.5), 50 mM NaCl, 0.1 mM zinc acetate and 0.01% Triton® X-100.

Unit Definition: One unit converts 1 µg of heat-denatured calf thymus DNA into an acid-soluble form in 1 minute at 37°C.¹¹

Activity Assay: The activity assay is performed in a 500 µl reaction containing 30 mM sodium acetate (pH 4.6), 50 mM NaCl, 1 mM zinc acetate, 500 µg/ml heat-denatured DNA and 0.01% Triton X-100.

10X Mung Bean Reaction Buffer: is 300 mM sodium acetate (pH 4.6), 500 mM NaCl, 10 mM zinc acetate and 0.1% Triton X-100.

Contaminating Activity Assays: Mung Bean Nuclease is free of detectable acid phosphatase and alkaline phosphatase activities. Double-stranded exo- and endonuclease activity is less than 0.05% of the enzyme's activity on single-stranded DNA.

Related Products: The following products are also available:

- Exonuclease I
- Exonuclease III
- Plasmid-Safe™ ATP-Dependent DNase
- RNase-Free DNase I
- T4 Polynucleotide Kinase
- EZ::TN™ Plasmid-Based Deletion Machine
- pWEB::TNC™ Deletion Cosmid Transposition Kit
- EZ::TN™ Transposon Insertion Kits
- SequiTherm™ Cycle Sequencing Kits

References:

1. Sung and Laskowski, M. (1962) *J. Biol. Chem.* **237**, 506.
2. Kowalski, D. *et al.*, (1976) *Biochemistry* **15**, 4457.
3. Sambrook, J. *et al.*, (1989) *Molecular Cloning: A Laboratory Manual (2nd ed.)*, New York: Cold Spring Harbor Laboratory Press.
4. Kroeker, W.D. *et al.*, (1976) *Biochemistry* **15**, 4463.
5. Green, M.R. and Roeder, R.G. (1980) *Cell* **22**, 231.
6. Mathis, D.J. *et al.*, (1981) *Proc. Natl. Acad. Sci. USA* **78**, 7383.
7. Murray, M.G. *et al.*, (1986) *Anal. Biochem.* **158**, 165.
8. Calzone, F.J. *et al.*, (1987) *Meth. Enzymol.* **152**, 611.
9. Henikoff, S. (1984) *Gene* **28**, 351.
10. Gubler, U. (1987) *Meth. Enzymol.* **152**, 330.
11. Vogt, V.M. (1973) *Eur. J. Biochem.* **33**, 192.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

Plasmid-Safe, EZ::TN, pWEB::TNC and SequiTherm are trademarks of EPICENTRE, Madison, Wisconsin.