

## TransforMax™ EC100™ Chemically Competent *E. coli*

Cat. No. CC02810

TransforMax™ EC100™ Chemically Competent *E. coli* are useful for primary cloning as well as subcloning applications.

The cells are provided in 50 µl aliquots (1 transformation per tube) for ease of use. The transformation efficiency of the TransforMax EC100 Chemically Competent cells is  $>1 \times 10^8$  cfu/µg DNA using pUC19.

### Relevant Phenotype

- Compatible with vectors expressing the LacZ' α-complementing peptide for "blue/white" screening of recombinants.
- Restriction minus for efficient cloning of methylated DNA.
- Accepts large clones.
- Endonuclease minus (*enda1*) to ensure high yields of plasmid clones.
- Recombination minus (*recA1*) to ensure the stability of large cloned inserts.

**Related Products:** The following products are also available:

- EZ-Tn5™ Transposons
- EZ-Tn5™ Transposome™ Kits
- TransforMax™ EC100™ Electrocompetent *E. coli*
- TransforMax™ EC100D™ *pir*<sup>+</sup> and *pir*-116 Electrocompetent *E. coli*
- Fast-Link™ DNA Ligation Kits

Cat. #

Quantity

---

**TransforMax™ EC100™ Chemically Competent *E. coli* includes:**

10 tubes of cells

50 µl/tube

pUC19 Control DNA (100 pg/µl)

10 µl

## Product Specifications

**Storage:** Store TransforMax EC100 *E. coli* cells at  $-70^{\circ}\text{C}$  and the pUC-19 Control DNA at either  $-20^{\circ}\text{C}$  or  $-70^{\circ}\text{C}$ . Do not thaw and refreeze the cells. Refreezing will result in significantly reduced transformation efficiency.

### Genotype:

F<sup>-</sup> *mcrA*  $\Delta$ (*mrr-hsdRMS-mcrBC*)  $\phi$ 80*dlacZ* $\Delta$ M15  $\Delta$ *lacX74* *recA1* *endA1* *araD139*  $\Delta$ (*ara, leu*)7697 *galU galK*  $\lambda$ -*rpsL nupG*.

**Quality Control:** TransforMax EC100 Chemically Competent *E. coli* yield  $>1 \times 10^8$  transformants per microgram of supercoiled pUC19 DNA. Transformation efficiency is determined using 10 pg of pUC19 and the Standard Transformation Procedure described on the back of this sheet. TransforMax EC100 Chemically Competent *E. coli* are tested to be free of contaminating DNA rendering resistance to ampicillin, tetracycline, kanamycin and chloramphenicol.

## Transformation Procedures

Two procedures for transforming the TransforMax EC100 Chemically Competent *E. coli* are presented. The Standard Transformation Procedure will provide the highest transformation efficiency. The 5 Minute Transformation Procedure is more rapid but may yield 10-fold or more lower transformation efficiency. The 5 Minute Transformation Procedure should only be used with ampicillin selection of the clones. Both procedures were written for transformation of 50  $\mu\text{l}$  of TransforMax EC100 Chemically Competent *E. coli*.

**Note:** Once thawed, do not refreeze the cells.

A different volume of cells can also be used based on the experiences and needs of the user.

## Standard Transformation Procedure

1. Prepare 250  $\mu\text{l}$  of SOC medium [Hanahan, D., (1983) *J. Mol. Biol.*, **166**, 557] for each transformation to be performed. Maintain the media at room temperature.
2. Heat a water bath or other temperature-controlled apparatus to  $42^{\circ}\text{C}$ .
3. Thaw the appropriate number of tubes of TransforMax EC100 Chemically Competent *E. coli* cells on ice. Mix by gentle tapping. Use the cells immediately.
4. Transfer 1-5  $\mu\text{l}$  of DNA or ligation reaction into each tube. Cap the tubes and incubate on ice for 5-30 minutes.
5. Transfer the tubes to  $42^{\circ}\text{C}$  and heat shock for 30 seconds.
6. Transfer the cells back to ice and cool for 2 minutes.
7. Remove the cover of the tubes and add 250  $\mu\text{l}$  of SOC Media.
8. Recover the cells by incubating at  $37^{\circ}\text{C}$  for 60 minutes with horizontal shaking (e.g. 225 rpm).
9. Plate the cells on the appropriate media and antibiotic, and grow overnight at  $37^{\circ}\text{C}$ .

### 5 Minute Transformation Procedure

The rapid 5 Minute Transformation Procedure may yield 10-fold or more lower transformation efficiency than the Standard Transformation Procedure described above. Importantly, only selection with ampicillin can be used with the 5 Minute Transformation Procedure.

1. Thaw the appropriate number of tubes of TransforMax EC100 Chemically Competent *E. coli* cells on ice. Mix by gentle tapping. Use the cells immediately.
2. Transfer 3-5  $\mu$ l of DNA or ligation reaction into each tube. Cap the tubes and incubate on ice for 5 minutes.
3. Spread the entire cell/DNA mixture onto a pre-warmed LB+ampicillin (100  $\mu$ g/ml) plate and grow overnight at 37°C.

*EC100, EC100D, EZ-Tn5, Transposome, Fast-Link, and TransforMax* are trademarks of Epicentre, Madison, Wisconsin.

Visit our technical blog: [epicentral.blogspot.com](http://epicentral.blogspot.com)

epicentre®  
an Illumina® company

726 Post Road, Madison, WI 53713 (800) 284-8474 (608) 258-3080 Fax (608) 258-3088